




ADVANCED PROCESS
MODELLING FORUM 2017
London 25–26 April


Introducing gPROMS FormulatedProducts

A comprehensive, integrated design tool
for formulated products manufacture

Dana Barrasso – Senior Consultant




A “new” PSE product gPROMS FormulatedProducts



gCRYSTAL, gSOLIDS, gCOAS

Adoption by industry



■ Active Ingredients


- Process development
- Manufacture


■ Formulations

- Product & process development
- Manufacture

■ Product performance


- Oral absorption
- *Product stability*
- *Dermal absorption*
- *Food digestion*
- *Agrochemical take-up by plants*
-






Tablet structure, e.g.

- Granule PSD
- Extra granular excipient
- Ungranulated API / excipient




Granule structure, e.g.

- PSD of primary particles (API & excipient)
- Porosity or liquid content




Primary particles, e.g.

- Impurities




gPROMS Formulated Products – Interoperability

Linking manufacture directly to product performance




**Mass of free API primary particles
suspended in tablet dissolver**





gPROMS Formulated Products – Interoperability

Linking manufacture directly to product performance




**Mass of free API primary particles
suspended in tablet dissolver**





Mass of free API primary particles suspended in tablet dissolver



gPROMS FormulatedProducts

Model libraries and features



Synthesis and Crystallization

Well-mixed reactor

Plug flow reactor

Well-mixed crystallizer

Plug flow crystallizer

Wet mill

Pressure filter

Filter dryer

Solids Processing

Dry mill

Roller compactor

Agglomerator

High shear wet granulator

Twin screw granulator

Fluid-bed granulation

2D pop balance granulation framework

Fluid-bed drying

Spray drying

Hot melt extruder

Lyophilization

Tablet press

Conveying

Tablet coater

Continuous blenders

Film coater

Bin flow / segregation

Product Performance

In-vitro dissolution and precipitation

Oral absorption in GI tract

Pharmacokinetics

Stability

Properties GUI

Data import and processing

Global Sensitivity Analysis

Parallelisation of GSA and PE

Multizonal models

gFP <-> DEM links

Extension of gSAFT databank

Transferred from gCRYSTAL, gSOLIDS or gCOAS

New in gPROMS FormulatedProducts MODEL In development

- Flexible database structure compatible with
 - PSE provided databases
 - 3rd party databases
 - corporate databases, etc.

Materials



Dosage forms



Equipment



Physiology



- Significant increase in usability
 - single repository for validated data
 - less looking up of data
 - fewer transcription errors

gPROMS Formulated Products – Data handling Database management



gPROMS Data Import & Physical Properties Configuration 1.0.0 - gPROMSDataImport

File Database Experiment Table Window Help

Edit Data import Properties database

PSE database

- Reactions
- Species
- Phase equilibria
 - Equilibrium adsorbed moisture content
 - Solubility
 - acetone_paracetamol_water
 - isopropanol_fenofibrate_water
 - l glutamic_acid_water
 - paracetamol_water
 - water_atenolol
 - water_fenofibrate
 - water_ketoprofen
 - water_lactose_solubility
 - water_plus_metoprolol
 - Solubility product
- Paracetamol synthesis and tableting
 - Phases
 - API crystals
 - Aqueous phase
 - Granules
 - Lactose excipient
 - Tablet
 - Phase transitions
 - API crystals:Aqueous phase

gPROMS Physical Properties package - (acetone_paracetamol_water)

Species in phase equilibrium

Solubility

Data entries

Solubility polynomial (anti-solvent)

Solubility polynomial

Model selected: Solubility polynomial (anti-solvent)

Solute species: paracetamol

Solvent species: acetone

Anti-solvent species: water

Mass fraction cutoff: 0

a: 3.682

b: 48.3866 1/(kg/kg liquid without solu...

c: 102.422 1/(kg/kg liquid without solu...

450.099 1/(kg/kg liquid without solu...

-576.746 1/(kg/kg liquid without solu...

-0.99207

-1401.0 K

Editing phase: API crystals

Phase name: API crystals

Phase type: Solid

Constituents: Single species

Constituent species: paracetamol

PhysProp package: gFPPP

Model: Basic model (solid)


Custom properties: No

OK Cancel

gPROMS Formulated Products – Data handling Data import and processing – previously



Pilotlab		Data Report																
Product name & ID		Product desc.						Batch Number			ID							
star 1.1		nr						1			784							
Recipe name		Report Filename						Batch started at			Batch ended at							
NIR V1		star 1.1_d.dat						04.11.14 11:36:20			04.11.14 13:01:13							
Step	Date/Time	TIR1.01	TIR3.01	TIR3.02	TIR3.03	TIR3.41	PICR3.0	PICR2.4	PICR2.2	PICR2.0	QIR2.01	MIR3.02	MIR3.41	QIR3.4	PDR3.41	QICR3.41		
	dd.mm.yy hh:mm	[°C]	[°C]	[°C]	[°C]	[°C]	[msh]	[bar]	[bar]	[g/min]	[g]	[%w]	[g/kg]	[%w]	[g/kg]	[%w]		
0	04.11.14 11:36:20	*** Process Start *** <NIR(Vers. 1)>																
1	04.11.14 11:36:20	Start "CONDITIONING" 1																
1	04.11.14 11:36:30	35.0	58.2	49.6	15.4	24.6	0	0.00	0.00	0.000	0.00	36.7	6.6	56.9	11.8	11.4	0.0	0.0
1	04.11.14 11:36:40	35.0	58.2	49.6	15.5	24.6	0	0.00	0.00	0.000	0.00	37.2	6.7	56.8	11.7	11.3	0.0	0.0
1	04.11.14 11:36:50	35.0	58.0	49.8	15.5	24.7	0	0.00	0.00	0.000	0.00	37.6	6.8	56.3	11.7	11.3	0.0	0.0
1	04.11.14 11:37:00	35.1	58.0	49.8	15.5	24.9	261	0.00	0.00	0.000	0.00	37.6	6.8	55.8	11.5	11.1	0.0	0.1
2	04.11.14 11:37:01	*** Process End *** <NIR(Vers. 1)>																
3	04.11.14 11:39:05	*** Process Start *** <NIR(Vers. 2)>																
4	04.11.14 11:39:06	Start "CONDITIONING" 1																
4	04.11.14 11:39:10	32.8	54.6	47.8	15.8	27.0	0	0.00	0.00	0.000	0.00	35.3	6.4	55.2	11.7	11.0	0.0	0.1
4	04.11.14 11:39:20	32.8	54.8	47.9	15.8	26.9	0	0.00	0.00	0.000	0.00	35.4	6.4	54.9	11.6	10.8	0.0	0.0
4	04.11.14 11:39:30	32.8	54.8	48.1	15.8	26.9	0	0.00	0.00	0.000	0.00	35.4	6.4	54.8	11.6	10.8	0.0	0.1
4	04.11.14 11:39:40	32.8	54.8	48.3	15.8	26.8	0	0.00	0.00	0.000	0.00	35.4	6.4	54.5	11.6	10.7	0.0	0.0
4	04.11.14 11:39:50	32.8	54.9	48.5	15.8	26.9	315	0.00	0.00	0.000	0.00	35.2	6.4	54.1	11.5	10.7	0.2	0.0
4	04.11.14 11:40:00	32.8	56.0	48.4	15.9	27.2	601	0.09	0.01	0.000	0.00	34.9	6.3	53.7	11.4	10.7	0.2	0.0
4	04.11.14 11:40:10	32.8	60.1	47.5	16.0	27.8	678	0.18	0.03	0.000	0.00	34.3	6.2	54.2	11.5	11.0	0.2	0.0
4	04.11.14 11:40:20	33.3	65.5	49.2	16.0	28.2	633	0.21	0.12	0.000	0.00	34.1	6.1	54.5	11.7	11.0	0.2	0.0
4	04.11.14 11:40:30	33.7	69.8	44.8	16.1	28.8	588	0.21	0.19	0.000	0.00	34.2	6.2	54.9	11.8	11.0	0.2	0.0
4	04.11.14 11:40:40	34.2	73.2	43.5	16.1	29.0	562	0.21	0.18	0.000	0.00	34.2	6.2	54.7	11.9	11.0	0.2	0.0
4	04.11.14 11:40:50	34.2	75.8	42.5	16.1	29.3	533	0.21	0.19	0.000	0.00	34.3	6.2	54.6	11.9	11.0	0.2	0.0
4	04.11.14 11:41:00	34.7	77.7	41.5	16.1	29.7	537	0.19	0.20	0.000	0.00	34.3	6.2	54.5	11.9	10.8	0.1	0.1
4	04.11.14 11:41:10	35.0	78.5	40.6	16.2	29.9	502	0.20	0.22	0.000	0.00	34.3	6.2	54.5	12.0	11.0	0.2	0.0
4	04.11.14 11:41:20	35.1	78.6	39.9	16.2	30.1	510	0.21	0.21	0.000	0.00	34.4	6.2	54.6	12.1	11.0	0.1	0.1
4	04.11.14 11:41:30	35.3	78.1	39.1	16.2	30.2	501	0.19	0.21	0.000	0.00	34.4	6.2	54.5	12.1	11.0	0.2	0.0
4	04.11.14 11:41:40	35.4	77.5	38.5	16.2	30.3	496	0.20	0.21	0.000	0.00	34.4	6.2	54.6	12.2	11.0	0.2	0.0
4	04.11.14 11:41:50	35.4	76.8	37.8	16.2	30.8	491	0.21	0.20	0.000	0.00	34.4	6.2	54.6	12.3	11.0	0.2	0.0
4	04.11.14 11:42:00	35.6	76.2	37.2	16.2	30.8	506	0.19	0.19	0.000	0.00	34.3	6.2	54.5	12.3	10.8	0.1	0.1
4	04.11.14 11:42:10	35.9	76.0	36.5	16.2	31.0	504	0.20	0.19	0.000	0.00	34.3	6.2	54.5	12.3	11.0	0.1	0.0
4	04.11.14 11:42:20	35.9	75.8	36.0	16.2	31.0	491	0.20	0.19	0.000	0.00	34.3	6.2	54.9	12.4	11.0	0.1	0.0



My_data (New project)

Time	Variable Name	Add new
	Sensor	
	Variance model	
Add new time		



Select measured variables... Delete row Delete column Transpose

General Controls Measured data gPROMS language Properties


- Creating performed experiment entities is a key step in the model validation workflow.
 - Significant manual effort has been required.
 - Often there are many experiments, and manual steps are repeated.
 - Data may not be in the right format, right units, same place.

gPROMS Formulated Products – Data handling

Data import and processing – now



1. Data file import and link to gPROMS variables



3. Export to gPROMS

Environmental inputs

- External disturbances
- Commercial environment


Innovate UK collaborative R&D project on
[Global Systems Analysis](#)
AstraZeneca, Britest, GSK, Pfizer, PSE
2014/06 – 2016/05

Variability


Decisions

- Design
- Operational

Model parameters




Uncertainty



KPIs

- CQAs
- Process operability
- Process safety
- Environmental impact
- Economic performance



Risk associated
with model-based
decisions

Transitioning to gPROMS Formulated Products



Is gPROMS FormulatedProducts 1.0 backward compatible with respect to gCRYSTAL 4.2, gSOLIDS 4.1, and gCOAS 1.3?

■ User specifications

- We currently only guarantee backward compatibility with respect to **physical properties** entered in via the global specifications objects.
- We are aiming for a future release to provide back compatibility with respect to specifications entered in **unit operations**.

■ Model results

- **Yes**
- PSE has performed extensive testing to ensure that the results obtained in gPROMS FormulatedProducts match those obtained in gCRYSTAL, gSOLIDS, and gCOAS.

When should I start using gPROMS Formulated Products?

- When creating a new flowsheet model
- When you'd like to use new functionality with an old flowsheet model
 - Adding a new unit operation that wasn't previously available or interfacing with the properties database.
 - Some manual conversion will be necessary.
 - Note that you will be able to use **Global System Analysis** and the **Data Import Tool** without migrating to the new libraries.


Will my licenses support gPROMS FormulatedProducts?

- **Yes** – Licenses for gCRYSTAL, gSOLIDS, and gCOAS will allow you to launch gPROMS FormulatedProducts and open and use the new versions of the libraries you previously could access.
- gPROMS FormulatedProducts will contain libraries for which you do not have licenses:
 - Libraries corresponding to previous products for which you did not have licenses.
 - New libraries containing models not existing in any previous product (reactor models, plug flow crystallizer, wet mill, continuous blenders, high shear and twin screw granulators)
- For the libraries which are new in gPROMS FormulatedProducts 1.0, PSE is happy to offer free 3-month evaluation licences to start no later than 1 June 2017.

In summary, gPROMS FormulatedProducts ...



- ... is PSE's new tool to deliver integrated design of formulated products and their manufacturing processes
- ... merges and builds on gCRYSTAL, gSOLIDS and gCOAS
- ... has been developed with input from
 - SbP Alliance: Pfizer, Lilly, PSE
 - Collaborative R&D projects, partially funded by UK government
 - D3P (Digital Design for Drug Products)
 - REMEDIES
 - ADDoPT



- Formulated Products team
- Software Technology Group
- Many others at PSE
- Formulated Products Advisory Board
- Collaborators across industry and academia
 - SbP Alliance
 - ADDoPT
 - REMEDIES
 - D3P

Thank you

