

Understanding Twin Screw Granulation using Global Systems Analysis


Gavin Reynolds, David Streater, Dan Davies

Advanced Process Modelling Forum, London, UK

April 2017






Introduction



[Seem *et al.* 2015, *Pow Tech* 276]



Introduction




- Highly configurable screw layout
- Main process parameters
 - Powder feed rate
 - Liquid feed rate
 - Screw speed

Finding an optimal design and operational space for a new formulation can consume significant time and material

Can modelling make this more efficient?




Introduction




- Many rate processes occur simultaneously during the granulation process
- Population balance modelling can be used to model the changing properties of the granules as a result of these rate processes
- gFormulate has a model for a twin screw granulator




gFormulate



TSG Compartmental Model Overview




TSG Model Input Parameters




Reducing the complexity – what do we know?

- Formulation properties
 - Density
 - Packing fraction
 - Particle Size Distribution
- Process parameters
 - Powder feed rate
 - Water feed rate
 - *Screw speed*




- Screw configuration




Reducing the complexity – what can we estimate?

- Simplify number of compartments (blocks rather than each element)




- Residence time distribution
 - Significantly longer in kneading sections
 - Estimated from published graphs
 - Kumar *et al.* 2014 *Eur. J. Pharm. Biopharm.* 87



[Seem *et al.* 2016, *Eur. J. Pharm. Biopharm.* 106]



Unknown rate parameters



Methodology to estimate rate parameters

- Utilise global systems analysis (GSA) to automatically interrogate model across parameter space
- Calculate a square error (SSE) of the granule size distribution with measured data


All parameters, large range

Identify feasible region for each parameter – Graphical Comparison and SSE


Eliminate less significant variables

5 parameters (3 breakage, 2 layering), small range

Find parameter sets that minimise SSE



Experimental Granule Size Distribution




- Increasing proportion of large granules when moving from low to high granulating water amount.




Eliminating Parameters – Nucleation

- **Drop size**

- Lower values produce a peak at low granule size that does not match experimental data
- Above 1.5 mm there is little to no effect



Eliminating parameters – Consolidation




$$\frac{d\varepsilon}{dt} = -k(\varepsilon - \varepsilon_{min})$$

Consolidation Rate Constant

Minimum porosity

- Consolidation rate**

- Very little effect
- This parameter primarily influences the granule porosity, rather than the granule size
- This parameter cannot be directly estimated from granule size data



Eliminating parameters – Layering



$$G = G_m \frac{M_{\text{powder}}}{k M_{\text{granulate}} + M_{\text{powder}}} \exp[-a(x_w - x_{wc})^2]$$

Max Growth Rate



Kinetic K

Kinetic A

Minimum Moisture Content

• Max Growth Rate

- Low values give an unusual bimodal distribution



• Kinetic K

- High values give an unusual bimodal distribution



Eliminating parameters – Layering



$$G = G_m \frac{M_{\text{powder}}}{k M_{\text{granulate}} + M_{\text{powder}}} \exp[-a(x_w - x_{wc})^2]$$



Max Growth Rate (points to G_m)

Kinetic K (points to k)

Kinetic A (points to a)

Minimum Moisture Content (points to x_{wc})


- **Kinetic A**
 - Little effect



- **Min moisture content**
 - Little effect except at very high values



Eliminating parameters – Breakage



Breakage Rate Constant



Alpha

$$S_M = \begin{cases} S c_M \left(\frac{y}{y_c} \right)^\alpha & \text{for } y \geq y_c \\ 0 & \text{for } y < y_c \end{cases}$$

Critical Particle Size

- **Breakage rate constant**

- Significant effect
- Low values give unrealistically large granules




- **Critical particle size**

- Unusual distribution above 500μm
- Below this considered physically unrealistic



Breakage rate and Growth rate – GSA ranked proximity to experimental data


Overall Breakage Rate



$$S_M = \begin{cases} S c_M \left(\frac{y}{y_c} \right)^\alpha & \text{for } y \geq y_c \\ 0 & \text{for } y < y_c \end{cases}$$

18


Overall Growth Rate



$$G = G_m \frac{M_{\text{powder}}}{k M_{\text{granulate}} + M_{\text{powder}}} \exp[-a(x_w - x_{wc})^2]$$



Breakage rate and experimental d50




Increasing Experimental D50

Decreasing Model Optimal Breakage Rate




Significance of breakage in TSG

Breakage in Distributive Mixing Elements




Breakage probability of material stronger than 9 kPa is independent of the dynamic yield strength



Daughter size distribution is independent of material strength



Tailored Granule Design



Different particle size using different screw design



'Investigating the effect of Raw Material Properties on the Breakage Rate Process in a Twin Screw Granulator',

S. Pradhan, M. Sen, C. Wassgren, J. Litster, I. Gabbott, G. Reynolds,
AIChE Annual Meeting 2016



Summary

- Understanding/Observations
 - Model complexity can be reduced by identifying parameters that don't have a significant effect on output properties
 - Some parameters can't be estimated by only using granule size data (e.g. consolidation rate)
 - Breakage appears to be a dominating mechanism in TSG for granule size
- Some areas for model development
 - Estimated rate parameters are a function of liquid amount
 - Ideally liquid amount would be mechanistically captured by the rate expressions
 - Screw speed is not explicitly captured – this could influence RTD in addition to collision frequency / shear stress





Confidentiality Notice

This file is private and may contain confidential and proprietary information. If you have received this file in error, please notify us and remove it from your system and note that you must not copy, distribute or take any action in reliance on it. Any unauthorized use or disclosure of the contents of this file is not permitted and may be unlawful. AstraZeneca PLC, 2 Kingdom Street, London, W2 6BD, UK, T: +44(0)20 7604 8000, F: +44 (0)20 7604 8151, www.astrazeneca.com

